



Professional Epilating System



2 Year Warranty

Precision Digital/Analog Technology

- Easy to read digital displays.
- Digital technology provides precision current application to assure
 - Increased client comfort.
 - Pure RF output, combines simplicity and effectiveness.
- One touch system displays DC, RF, and timer settings.
- Two footswitches or one footswitch for blend operation.
- Manual adjustments allows electrologist full control with single current (RF or DC) or blend modality.
- Features pre-programmed selections for various hair types however settings are completely manually adjustable by electrologist.
- Phoresis: CATA+ post treatment for electrolysis. ANA- for skin care.
- Timer: Audio (soft tone) for settings and treatment timer (1-20 seconds).
- Internal power supply is highly energy efficient and adapts itself to all worldwide AC power standards.

Clareblend Inc. is an ISO13485 certified and FDA registered manufacturer. All Clareblend equipment is designed and manufactured in USA at our Reno, Nevada facility.



Clareblend 3555 Airway Dr. #307 Reno NV 89511 Toll Free: 800-334-7126 Local: 775-332-3850 Email: staff@clareblend.com

