



3 Year Warranty

Newest and latest technology to provide faster response, extreme accuracy, and precise timing. This assures maximum client comfort.

The Elegance is a <u>true blend</u> epilator

- Features auto memory recall of previously used settings when switching modalities. Four programmable custom settings to store favorite treatment settings.
- Proven quick release presets for all hair types.
 Preset settings remain fully adjustable by the electrologist.
- Air feature available in Blend & single currents modalities soothes the surface of the skin.
 Selectable by foot switch activation, autosense or set for continuous flow.
- Highly efficient internal power supply saves energy for low cost operation.

Operating Features:

Single and BLEND current modalities:

RF, DC or full BLEND:

These currents can be activated by:

- Two footswitches for blend or single current (RF or DC).
- One footswitch (Autoblend mode) for blend or single current (RF or DC).
- No footswitch (Autosense mode) for single or dual currents.



Additional Elegance features:

Pulse mode:

 RF pulsing. Available with blend or RF modalities for additional client comfort.

UL display:

 Units of lye automatically calculated when using DC current and treatment timer.

Timers:

• Two timers for treatment (0-20 seconds) and session (5-60 minutes) timing.

Phoresis:

- CATA+ for post electrology treatments.
- ANA- for skin care treatments.

Audio:

 Audio feedback (soft tone) for operator settings and timing functions.

Accessory testing:

 Test feature for testing probe cord, footswitches and accessories.

The Elegance is always in your control!

Clareblend Inc. is an ISO13485 certified and FDA registered manufacturer.

All Clareblend equipment is designed and manufactured in USA at our Reno, Nevada facility.





Clareblend Inc. 3555 Airway Dr. #307 Reno NV 89511 Toll Free: 800-334-7126 Local: 775-332-3850 Email: staff@clareblend.com